

Dealers Stamp

Dan Medica South Ltd
Exclusive Supplier of

DebbonAir
Dynamics

High Risk
Wheel Chair Cushions
& Mattresses

Dan Medica South Ltd

28 Downsview Avenue, Storrington, West Sussex RH20 4PS

Tel: **00 44 208 133 2851**

Fax: **00 44 208 196 2364**

e-mail: info@danmedicasouth.co.uk

web: www.danmedicasouth.co.uk

Warranty


Each and every DebbonAir Dynamics cushion and mattress is individually serialised. This allows us and your supplier to trace the exact manufacturing and supply date along with detailed information held in our factory with respect to materials, etc. Your DebbonAir Dynamics product has a manufacturers limited guarantee of 2 years from date of purchase

DanMedica South Limited are proud to be the exclusive supplier of the **DebbonAir Dynamics** Range of cushions and mattresses. manufactured to combat even the highest risk of pressure sores.

History and Manufacturing

Over the last 6 years many different materials have been considered and indeed tried, ranging from the more common neoprene and rubber used in products that have been marketed for many years, these proved to have no benefit in our manufacturing process which is very different over and above our competitors. Instead of antiquated dipping with neoprene or rubber, our **DebbonAir** products have NO separate airway channels, each product is Injection moulded plus we use High Frequency Welding (patented) to bring back and front sections together. There is NO use of glue or any other bonding agents applied to our products at manufacturing, hence we do not suffer back separations of any kind or the related problems used when glue is required

Our very special PVC compounds used in both top and bottom sections are of course fire/flamm retardant, anti allergic to direct skin contact, non fungal and environment friendly. Our base section PVC is so strong it is currently used in sea conditions to repel oil spillage.



The DebbonAir cell walls have been designed far deeper than any of our competitors, our products are all blue in colour to distinguish **DebbonAir** over and above any other manufacturers who market finned cushions.

Our patented High Frequency welding using the most up to date PVC materials brings the DebbonAir product range into the market at a far more competitive price, due to this latest technology.

Available Sizes

The DebbonAir 6cm height cushion is available in 6 standard sizes with no price difference between one or two air valves. The DebbonAir 10cm is also available in 6 standard sizes, also without any price difference for one or two air valves.

DebbonAir Cushions

	6cm	10cm	1 valve	2 valve
36cm x 36cm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
36cm x 40cm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
40cm x 40cm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
44cm x 40cm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
48cm x 44cm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
48cm x 48cm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



Cushion Covers

Please note all DebbonAir Dynamic products, excluding mattresses, have full two way stretch fire retardant covers, and are fully compliant with British and European regulations (CE mark).



Special Sizes

Due to the most advanced manufacturing techniques in the world, we are able to supply almost any size cushion or mattress. We are also able to fit DebbonAir products into special seating, moulded seats, wheelchair backs, lift & recliner chairs plus a host of other pressure care equipment where only the latest technology will suffice. CPR valves can be fitted to all our products at an extra charge.


We are pleased to inform both our direct sales customers and our distributors, costs involved when manufacturing specific custom made products shall be kept to a minimum in utilising this innovative new technology Current delivery on special products is approximately 2 weeks.

DebbonAirMaitress Sections

(note 4 sections per standard bed)

	6cm	10cm	1 valve	2 valve
100cm x 50cm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

The above sizes are available in one or two air valves. When ordering please state size, cell height and number of air valves.



Advantages of the New DebbonAir Dynamics product range

- Available in any size
- No glue or other bonding agents used
- Patented high frequency welding
- Manufactured from tough PVC materials
- The very latest Injection Moulding Technology
- No high or low spots in manufacturing
- Less likely to puncture
- Same outlay for 1 or 2 air valves
- Same outlay for 6cm or 10cm heights
- Full 2 Year manufacturers warranty
- Blue finned deep air cell construction
- New integrated Fast air flow system
- C.P.R exhaust valves available on all products
- Two way airpumps for fast inflation
- Complete with cover, repair kit, pump & manual
- Supplied in presentation box

Dan Medica South Ltd



28 Downsview Avenue
Storrington,
West Sussex, RH20 4PS

Tel: 00 44 208 133 2851
Fax: 00 44 208 196 2364

email: info@danmedicasouth.co.uk
web: www.danmedicasouth.co.uk

UK Reg No. 437 3372